

Product datasheet for **GA107831**

Human PHLDA1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHLDA1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology like domain family A member 1
Symbol:	PHLDA1
Components:	GA107831G1, PHLDA1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107831G2, PHLDA1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107831G3, PHLDA1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q8WV24
Synonyms:	DT1P1B11; PHRIP; TDAG51
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)