

Product datasheet for **GA107872**

Human PEPP3 (PLEKHA6) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLEKHA6 CRISPRa kit - CRISPR gene activation of human pleckstrin homology domain containing A6
Symbol:	PLEKHA6
Components:	GA107872G1, PEPP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107872G2, PEPP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107872G3, PEPP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9Y2H5
Synonyms:	PEPP-3; PEPP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)