

Product datasheet for **GA107913**

Human KIFAP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KIFAP3 CRISPRa kit - CRISPR gene activation of human kinesin associated protein 3
Symbol:	KIFAP3
UniProt ID:	Q92845
Synonyms:	dj190I16.1; FLA3; KAP-1; KAP-3; KAP3; SMAP; Smg-GDS
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)