

Product datasheet for **GA107923**

Human POMZP3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	POMZP3 CRISPRa kit – CRISPR gene activation of human POM121 and ZP3 fusion
Symbol:	POMZP3
Components:	GA107923G1, POMZP3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107923G2, POMZP3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107923G3, POMZP3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PJE2</u>
Synonyms:	POM-ZP3; POM121
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

