


Product datasheet for **GA107934**

Human Kanadaptin (SLC4A1AP) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SLC4A1AP CRISPRa kit - CRISPR gene activation of human solute carrier family 4 member 1 adaptor protein
Symbol:	SLC4A1AP
Components:	GA107934G1, Kanadaptin gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107934G2, Kanadaptin gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107934G3, Kanadaptin gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9BWU0</u>
Synonyms:	HLC3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)