

## Product datasheet for **GA107935**

### Human P2X2 (P2RX2) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	P2RX2 CRISPRa kit - CRISPR gene activation of human purinergic receptor P2X 2
Symbol:	P2RX2
UniProt ID:	<u><a href="#">Q9UBL9</a></u>
Synonyms:	DFNA41; P2X2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)