

Product datasheet for **GA107937**

Human **SCMH1** activation kit by **CRISPRa**

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | SCMH1 CRISPRa kit - CRISPR gene activation of human Scm polycomb group protein homolog 1 |
| Symbol: | SCMH1 |
| UniProt ID: | <u>Q96GD3</u> |
| Synonyms: | Scml3 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)