

Product datasheet for **GA107939**

Human PTIP (PAXIP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PAXIP1 CRISPRa kit - CRISPR gene activation of human PAX interacting protein 1
Symbol:	PAXIP1
Components:	GA107939G1, PTIP gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107939G2, PTIP gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107939G3, PTIP gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6ZW49</u>
Synonyms:	CAGF28; CAGF29; PACIP1; PAXIP1L; PTIP; TNRC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)