

Product datasheet for **GA107967**

Human **STK38L** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	STK38L CRISPRa kit - CRISPR gene activation of human serine/threonine kinase 38 like
Symbol:	STK38L
UniProt ID:	<u>Q9Y2H1</u>
Synonyms:	NDR2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)