

Product datasheet for **GA107977**

Human SEMCAP3 (PDZRN3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDZRN3 CRISPRa kit - CRISPR gene activation of human PDZ domain containing ring finger 3
Symbol:	PDZRN3
UniProt ID:	Q9UPQ7
Synonyms:	LNX3; SEMACAP3; SEMCAP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)