


Product datasheet for **GA107985**

Human DOP1A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DOP1A CRISPRa kit - CRISPR gene activation of human DOP1 leucine zipper like protein A
Symbol:	DOP1A
Components:	GA107985G1, DOP1A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA107985G2, DOP1A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA107985G3, DOP1A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q5JWR5
Synonyms:	dj202D23.2; KIAA1117
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)