

Product datasheet for **GA107999**

Human SMG1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SMG1 CRISPRa kit - CRISPR gene activation of human SMG1 nonsense mediated mRNA decay associated PI3K related kinase
Symbol:	SMG1
UniProt ID:	Q96Q15
Synonyms:	61E3.4; ATX; KIAA0421; LIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)