

Product datasheet for **GA108015**

Human TXNDC4 (ERP44) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ERP44 CRISPRa kit - CRISPR gene activation of human endoplasmic reticulum protein 44
Symbol:	ERP44
UniProt ID:	Q9BS26
Synonyms:	PDIA10; TXNDC4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)