

Product datasheet for **GA108017**

Human KIAA0701 (UHRF1BP1L) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UHRF1BP1L CRISPRa kit - CRISPR gene activation of human UHRF1 binding protein 1 like
Symbol:	UHRF1BP1L
Components:	GA108017G1, KIAA0701 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108017G2, KIAA0701 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108017G3, KIAA0701 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A0JNW5
Synonyms:	SHIP164
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)