

Product datasheet for **GA108018**

Human SWAP70 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	SWAP70 CRISPRa kit - CRISPR gene activation of human switching B cell complex subunit SWAP70
Symbol:	SWAP70
UniProt ID:	Q9UH65
Synonyms:	HSPC321; SWAP-70
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)