

Product datasheet for **GA108034**

Human KIF1B activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | KIF1B CRISPRa kit - CRISPR gene activation of human kinesin family member 1B |
| Symbol: | KIF1B |
| Components: | GA108034G1, KIF1B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108034G2, KIF1B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108034G3, KIF1B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID: | <u>O60333</u> |
| Synonyms: | CMT2; CMT2A; CMT2A1; HMSNII; KLP; NBLST1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online »](#)