

Product datasheet for **GA108046**

Human TNRC6B activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNRC6B CRISPRa kit - CRISPR gene activation of human trinucleotide repeat containing adaptor 6B
Symbol:	TNRC6B
Components:	GA108046G1, TNRC6B gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108046G2, TNRC6B gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108046G3, TNRC6B gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9UPQ9
Synonyms:	KIAA1093
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)