

Product datasheet for **GA108049**

Human **TOGARAM1** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	TOGARAM1 CRISPRa kit – CRISPR gene activation of human TOG array regulator of axonemal microtubules 1
Symbol:	TOGARAM1
Components:	GA108049G1, TOGARAM1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108049G2, TOGARAM1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108049G3, TOGARAM1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y4F4</u>
Synonyms:	KIAA0423
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

