

Product datasheet for **GA108058**

Human Plexin D1 (PLXND1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PLXND1 CRISPRa kit - CRISPR gene activation of human plexin D1
Symbol:	PLXND1
UniProt ID:	Q9Y4D7
Synonyms:	PLEXD1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)