

Product datasheet for **GA108070**

Human DCUN1D4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DCUN1D4 CRISPRa kit - CRISPR gene activation of human defective in cullin neddylation 1 domain containing 4
Symbol:	DCUN1D4
Components:	GA108070G1, DCUN1D4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108070G2, DCUN1D4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108070G3, DCUN1D4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q92564
Synonyms:	FLJ42355; KIAA0276
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)