

Product datasheet for **GA108082**

Human WDR43 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WDR43 CRISPRa kit - CRISPR gene activation of human WD repeat domain 43
Symbol:	WDR43
Components:	GA108082G1, WDR43 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108082G2, WDR43 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108082G3, WDR43 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q15061</u>
Synonyms:	NET12; UTP5
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)