

Product datasheet for **GA108086**

Human MPRIP activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	MPRIP CRISPRa kit - CRISPR gene activation of human myosin phosphatase Rho interacting protein
Symbol:	MPRIP
UniProt ID:	<u>Q6WQC1</u>
Synonyms:	M-RIP; MRIP; p116Rip; RHOIP3; RIP3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)