

Product datasheet for **GA108089**

Human EFR3A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	EFR3A CRISPRa kit – CRISPR gene activation of human EFR3 homolog A
Symbol:	EFR3A
Components:	GA108089G1, EFR3A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108089G2, EFR3A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108089G3, EFR3A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q14156</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

