


## Product datasheet for **GA108091**

### Human **SLC35D1** activation kit by **CRISPRa**

#### Product data:

Product Type:	CRISPRa Kits
Description:	SLC35D1 CRISPRa kit - CRISPR gene activation of human solute carrier family 35 member D1
Symbol:	SLC35D1
Components:	GA108091G1, SLC35D1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108091G2, SLC35D1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108091G3, SLC35D1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q9NTN3</a></u>
Synonyms:	UGTREL7
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)