

Product datasheet for **GA108095**

Human Methionyl Aminopeptidase 1 (METAP1) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	METAP1 CRISPRa kit - CRISPR gene activation of human methionyl aminopeptidase 1
Symbol:	METAP1
Components:	GA108095G1, Methionyl Aminopeptidase 1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108095G2, Methionyl Aminopeptidase 1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108095G3, Methionyl Aminopeptidase 1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	P53582
Synonyms:	MAP1A; MetAP1A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)