

## Product datasheet for **GA108100**

### Human PASK activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PASK CRISPRa kit - CRISPR gene activation of human PAS domain containing serine/threonine kinase
Symbol:	PASK
UniProt ID:	<a href="#">Q96RG2</a>
Synonyms:	PASKIN; STK37
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)