

Product datasheet for **GA108107**

Human PHLDB1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PHLDB1 CRISPRa kit - CRISPR gene activation of human pleckstrin homology like domain family B member 1
Symbol:	PHLDB1
Components:	GA108107G1, PHLDB1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108107G2, PHLDB1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108107G3, PHLDB1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q86UU1</u>
Synonyms:	LL5A
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)