

Product datasheet for **GA108117**

Human Proteasome Activator Subunit 4 (PSME4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PSME4 CRISPRa kit - CRISPR gene activation of human proteasome activator subunit 4
Symbol:	PSME4
UniProt ID:	Q14997
Synonyms:	PA200
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)