

Product datasheet for **GA108120**

Human FAM168A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FAM168A CRISPRa kit - CRISPR gene activation of human family with sequence similarity 168 member A
Symbol:	FAM168A
Components:	GA108120G1, FAM168A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108120G2, FAM168A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108120G3, FAM168A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q92567</u>
Synonyms:	KIAA0280; TCRP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)