

## Product datasheet for **GA108145**

### Human SEL1L3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	SEL1L3 CRISPRa kit - CRISPR gene activation of human SEL1L family member 3
Symbol:	SEL1L3
UniProt ID:	<a href="#">Q68CR1</a>
Synonyms:	Sel-1L3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)