

## Product datasheet for **GA108150**

### Human Phospholipase C beta 1 (PLCB1) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PLCB1 CRISPRa kit - CRISPR gene activation of human phospholipase C beta 1
Symbol:	PLCB1
UniProt ID:	<a href="#">Q9NQ66</a>
Synonyms:	EIEE12; PI-PLC; PLC-154; PLC-I; PLC154; PLCB1A; PLCB1B
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)