

Product datasheet for **GA108153**

Human TMEM131L (KIAA0922) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TMEM131L CRISPRa kit - CRISPR gene activation of human transmembrane 131 like
Symbol:	TMEM131L
Components:	GA108153G1, TMEM131L gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108153G2, TMEM131L gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108153G3, TMEM131L gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2VDJ0</u>
Synonyms:	TMEM131L
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)