

Product datasheet for **GA108157**

Human SCC112 (PDS5A) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PDS5A CRISPRa kit - CRISPR gene activation of human PDS5 cohesin associated factor A
Symbol:	PDS5A
Components:	GA108157G1, SCC112 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108157G2, SCC112 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108157G3, SCC112 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q29RF7</u>
Synonyms:	PIG54; SCC-112; SCC112
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)