

Product datasheet for **GA108164**

Human OTUD3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	OTUD3 CRISPRa kit – CRISPR gene activation of human OTU deubiquitinase 3
Symbol:	OTUD3
Components:	GA108164G1, OTUD3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108164G2, OTUD3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108164G3, OTUD3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T2D3</u>
Synonyms:	DUBA4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

