

Product datasheet for **GA108169**

Human DENND5A activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DENND5A CRISPRa kit - CRISPR gene activation of human DENN domain containing 5A
Symbol:	DENND5A
Components:	GA108169G1, DENND5A gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108169G2, DENND5A gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108169G3, DENND5A gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6IQ26
Synonyms:	RAB6IP1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)