

Product datasheet for **GA108170**

Human DDHD2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	DDHD2 CRISPRa kit - CRISPR gene activation of human DDHD domain containing 2
Symbol:	DDHD2
Components:	GA108170G1, DDHD2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108170G2, DDHD2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108170G3, DDHD2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O94830</u>
Synonyms:	iPLA(1)gamma; SAMWD1; SPG54
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)