

Product datasheet for **GA108172**

Human PPIP5K2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	PPIP5K2 CRISPRa kit - CRISPR gene activation of human diphosphoinositol pentakisphosphate kinase 2
Symbol:	PPIP5K2
Components:	GA108172G1, PPIP5K2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108172G2, PPIP5K2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108172G3, PPIP5K2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q43314
Synonyms:	CFAP160; HISPPD1; IP7K2; VIP2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)