

Product datasheet for **GA108179**

Human TSPYL4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TSPYL4 CRISPRa kit - CRISPR gene activation of human TSPY like 4
Symbol:	TSPYL4
Components:	GA108179G1, TSPYL4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108179G2, TSPYL4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108179G3, TSPYL4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UJ04</u>
Synonyms:	dj486l3.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)