

## Product datasheet for **GA108180**

### Human CAMSAP1L1 (CAMSAP2) activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                        |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                           |
| Description:  | CAMSAP2 CRISPRa kit - CRISPR gene activation of human calmodulin regulated spectrin associated protein family member 2                                                                                                                                                                                                 |
| Symbol:       | CAMSAP2                                                                                                                                                                                                                                                                                                                |
| Components:   | GA108180G1, CAMSAP1L1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108180G2, CAMSAP1L1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA108180G3, CAMSAP1L1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q08AD1</a>                                                                                                                                                                                                                                                                                                 |
| Synonyms:     | CAMSAP1L1                                                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                                |

#### Product images:



[View online »](#)