

Product datasheet for **GA108181**

Human TASOR activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TASOR CRISPRa kit - CRISPR gene activation of human transcription activation suppressor
Symbol:	TASOR
Components:	GA108181G1, TASOR gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108181G2, TASOR gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108181G3, TASOR gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9UK61</u>
Synonyms:	C3orf63; RAP140; se89-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)