

Product datasheet for **GA108189**

Human LPHN3 (ADGRL3) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	ADGRL3 CRISPRa kit - CRISPR gene activation of human adhesion G protein-coupled receptor L3
Symbol:	ADGRL3
UniProt ID:	Q9HAR2
Synonyms:	CIRL3; CL3; LEC3; LPHN3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)