

Product datasheet for **GA108198**

Human BICD2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | BICD2 CRISPRa kit - CRISPR gene activation of human BICD cargo adaptor 2 |
| Symbol: | BICD2 |
| UniProt ID: | <u>Q8TD16</u> |
| Synonyms: | bA526D8.1; SMALED2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)