


Product datasheet for **GA108206**

Human FKBP15 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	FKBP15 CRISPRa kit - CRISPR gene activation of human FKBP prolyl isomerase 15
Symbol:	FKBP15
Components:	GA108206G1, FKBP15 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108206G2, FKBP15 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108206G3, FKBP15 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q5T1M5</u>
Synonyms:	FKBP133; KIAA0674; PPP1R76
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)