

## Product datasheet for **GA108209**

### Human hCAP D3 (NCAPD3) activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	NCAPD3 CRISPRa kit - CRISPR gene activation of human non-SMC condensin II complex subunit D3
Symbol:	NCAPD3
UniProt ID:	<a href="#">P42695</a>
Synonyms:	CAP-D3; CAPD3; hCAP-D3; hcp-6; hHCP-6
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)