


Product datasheet for **GA108220**

Human WASH complex subunit 7 (WASHC4) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	WASHC4 CRISPRa kit - CRISPR gene activation of human WASH complex subunit 4
Symbol:	WASHC4
Components:	GA108220G1, WASH complex subunit 7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108220G2, WASH complex subunit 7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108220G3, WASH complex subunit 7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q2M389
Synonyms:	MRT43; SWIP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)