

Product datasheet for **GA108224**

Human TBCID30 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TBCID30 CRISPRa kit – CRISPR gene activation of human TBC1 domain family member 30
Symbol:	TBCID30
Components:	GA108224G1, TBCID30 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108224G2, TBCID30 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108224G3, TBCID30 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9Y2I9</u>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

