


Product datasheet for **GA108240**

Human KIAA0323 (KHNYN) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	KHNYN CRISPRa kit - CRISPR gene activation of human KH and NYN domain containing
Symbol:	KHNYN
Components:	GA108240G1, KIAA0323 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA108240G2, KIAA0323 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA108240G3, KIAA0323 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O15037</u>
Synonyms:	KIAA0323
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)