


Product datasheet for **GA108241**

Human UBR4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	UBR4 CRISPRa kit - CRISPR gene activation of human ubiquitin protein ligase E3 component n-recognin 4
Symbol:	UBR4
UniProt ID:	Q5T4S7
Synonyms:	p600; RBAF600; ZUBR1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)