

## Product datasheet for **GA108250**

### Human PSD3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	PSD3 CRISPRa kit - CRISPR gene activation of human pleckstrin and Sec7 domain containing 3
Symbol:	PSD3
UniProt ID:	<u><a href="#">Q9NYI0</a></u>
Synonyms:	EFA6D; EFA6R; HCA67
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)