

Product datasheet for **GA108258**

Human TENC1 (TNS2) activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	TNS2 CRISPRa kit - CRISPR gene activation of human tensin 2
Symbol:	TNS2
UniProt ID:	Q63HR2
Synonyms:	C1-TEN; C1TEN; TENC1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)